

Light Sources

L9000 LED Light Source

- » LED technology
- » No need to replace bulbs
- » Designed to produce less environmental waste
- » Safelight cable compatible



L9000 Light Source

220-210-000i L9000 LED Light Source

6000-001-020 Firewire Cable (*not included with L9000, but needed to control light source from 1288 camera head*)



X8000 Light Source

220-200-000E X8000 Xenon Light Source

Light Source Accessories

Safelight Fiberoptic Light Cables and Adapters



- 233-050-199** Safelight Kit:
Safelight Cable (233-050-100) & Safelight Scope Adapter (233-050-087)
- 233-050-100** 5.0 mm x 10ft Clear-Case Safelight Fiberoptic Light Cable
- 233-050-200** 5.0 mm x 10ft Opaque-Case Safelight Fiberoptic Light Cable
- 233-050-087** Safelight Henke/Dyonics/Storz Scope Adapter
- 233-050-088** Safelight Wolf Scope Adapter
- 233-050-086** Safelight Olympus Scope Adapter
- 233-050-192** Safelight Circon Scope Adapter



233-065-010 233-050-090 233-050-064 233-050-065 233-050-069 233-050-067

Fiberoptic Light Cables

- 233-065-010** 6.5 mm x 10ft Fiberoptic Light Cable
(Wolf Scope & Stryker Light Source Fittings)
- 233-050-090** 5.0 mm x 10ft Fiberoptic Light Cable
(Universal Scope and Stryker Light Source Fittings)
- 233-050-064** 5.0 mm x 10ft Clear Case Fiberoptic Light Cable
(Storz Scope & Stryker Light Source Fittings, Screw-On)
- 233-050-065** 5.0 mm x 10ft Fiberoptic Light Cable
(Stryker Scope & Stryker Light Source Fittings)
- 233-050-069** 5.0 mm x 10ft Clear Case Fiberoptic Light Cable
(Wolf Scope & Stryker Light Source Fittings, Snap-On)
- 233-050-067** 2.0 mm x 7.5ft Clear Case Fiberoptic Light Cable
(Stryker Scope & Stryker Light Source Fittings)
- 233-050-084** 5.0 mm x 10ft Non-Transparent Fiberoptic Cable

Replacement Bulbs



- 220-201-000** Replacement Bulb for X8000 (300 Watts)
- 220-191-000** Replacement Bulb for X7000 (300 Watts)

Light Source Adapters

Connect to Light Source
(For use with universal fiberoptic
light cable, 233-050-090)

- 233-050-071** Stryker/ACMI/Zimmer Light Source Adapter
- 233-050-072** Olympus Light Source Adapter
- 233-050-073** Storz Light Source Adapter
- 233-050-074** Wolf/Dyonics Light Source Adapter



233-050-075



233-050-076



233-050-077



233-060-000

Scope Adapters

Connect to Scope
(For use with universal fiberoptic
light cable, 233-050-090)

- 233-050-075** Stryker/ACMI Scope Adapter
- 233-050-076** Wolf/Dyonics Scope Adapter
- 233-050-077** Storz/Olympus Scope Adapter
- 233-060-000** ACMI Universal Scope Adapter



233-080-000



233-081-000



233-082-000

**Scope Adapters for
Non-Autoclavable Scopes**
(Black Eyepiece)

Connect to Scope

- 233-080-000** Stryker Adapter
- 233-081-000** Wolf/Dyonics Adapter
- 233-082-000** Storz Adapter



331-000-100



331-000-200



331-000-300



331-000-400

**Scope Adapters for
Autoclavable Scopes**
(Blue & Red Eyepiece)

Connect to Scope

- 331-000-100** Stryker Adapter
- 331-000-200** Wolf/Dyonics Adapter
- 331-000-300** Storz Adapter
- 331-000-400** Right Angle Adapter

*This does not replace the operative technique, please refer to the operative technique for all procedures. **Not drawn to scale.